```
In [1]:
            1 import warnings
           2 from pyIClab import (
                  DSM_CEConstrutor, IonChromatograph, Eluent, SwitchingValve,
            4
                  Column, IonExchanger, SampleLoop, Detector, Dummy,
            5
            6 | from pyIClab import ContaminatedPhreeqcSuppressorBeta as Suppressor
           7 from pyIClab.engines.equilibriums import find_x_LSSM
           8 from pyIClab.beadedbag import mpl_custom_rcconfig
           9 import matplotlib.pyplot as plt
           10 import seaborn as sns
           11 import numpy as np
           12 from IPython.display import clear_output
           13
           14 sns.set()
           15 | plt.rcParams.update(mpl_custom_rcconfig)
In [2]: ▼
           1 class LocalConstructor(DSM_CEConstrutor):
                  def set_x(self):
            3
            4
            5
                      kmap = self.set_kmap()
            6
            7
                      return find_x_LSSM(kmap, -1)
            1 | Vinj = ['25 uL', '100 uL', '500 uL', '1000 uL', '2000 uL']
In [3]:
In [4]:
            1 ic_collection = []
              for i, V in enumerate(Vinj):
            3
                  eluent = Eluent.HydroxideIsocratic('18.75 mM', name=f'EG{i}')
            4
                  sp = IonExchanger.load('home_made.dat', directory='db')
                  column = Column(f'Home_made{i}', length='15 cm', ID='4.6 mm')
            5
                  column.pack(sp)
            6
            7
                  sixport = SwitchingValve.SixPort(name=f'Sixport{i}')
            8
                  suppressor = Suppressor(name=f'Suppressor{i}', kind='anion', _CO2_level=.018)
            9
                  detector = Detector(name=f'Detector{i}')
           10
                  loop = SampleLoop(name=f'Loop{i}', V=V)
           11
           12
                  sixport.assemble(0, eluent)
           13
                  sixport.assemble([2, 5], loop)
           14
                  sixport.assemble(1, column)
           15
                  column.assemble(suppressor)
           16
                   suppressor.assemble(detector)
           17
           18
                  ic = IonChromatograph(name=f'IC{i}', competing_ions=('OH[-1]',), lockon=sixport)
           19
                  ic_collection.append(ic)
           20
           21 ic_collection
Out[4]: [<IC System "ICO">,
         <IC System "IC1">,
         <IC System "IC2">,
         <IC System "IC3">,
         <IC System "IC4">]
In [5]: ▼
           1 | for i, ic in enumerate(ic_collection):
                  solution = {'Cl-': '0.1 mM'}
            3
            4
                  ic.inject(solution, f'loop{i}')
                  commands = f''
            6
                  0.0 min, sixport{i}, inject
            7
            8
                  ic.reset_commands(commands)
                  ic.set_ModelConstructor(LocalConstructor, f'Home_made{i}')
            9
                  with warnings.catch_warnings(action='ignore'):
           10
                      ic.start(tmax='10 min')
           11
                      clear_output()
           12
```

```
In [6]:
            1 df_collection = []
            2 | for i, ic in enumerate(ic_collection):
            3
            4
                   (detector,) = tuple(ic.detectors)
            5
                   df_collection.append(detector.get_signals(signal_type='conductivity'))
             6
         Calculating eluent conductivity on <Detector "Detector0">...:
                                                                          0%
                                                                                       | 0/5999 [00:00<?, ?it/s]
         Calculating eluent conductivity on <Detector "Detector1">...:
                                                                          0%|
                                                                                       | 0/5999 [00:00<?, ?it/s]
         Calculating eluent conductivity on <Detector "Detector2">...:
                                                                          0%
                                                                                       | 0/5999 [00:00<?, ?it/s]
                                                                          0%|
                                                                                       | 0/5999 [00:00<?, ?it/s]
         Calculating eluent conductivity on <Detector "Detector3">...:
                                                                                       | 0/5999 [00:00<?, ?it/s]
         Calculating eluent conductivity on <Detector "Detector4">...:
                                                                          0%|
In [7]:
            1 fig, ax = plt.subplots()
             2 offsets = [0, 0, 0, 0, 0, 0]
            3 for i, df in enumerate(df_collection):
                   ax.plot(df['time'], df['signal']+offsets[i])
            5
                     ax.set(ylim=(0, 1))
            6
           400
           300
           200
           100
              0
                     0
                                       2
                                                                                             8
                                                         4
                                                                           6
                                                                                                              10
In [8]:
            1 | from scipy.signal import find_peaks
In [9]:
            1 |tR = []
            2
              for df in df_collection:
            3
                   (pk, *_), _ = find_peaks(df['signal']-df['signal'][0], height=2)
            4
                   tR.append(df['time'][pk])
            5
             6 tR
Out[9]: [4.951155385128376,
          5.0428434478159385,
          5.516287262420807,
          6.08975441813938,
          7.196679393131044]
            for i, df in enumerate(df_collection):
In [10]: ▼
            2
            3
                   df.to_csv(f'effective charge ratio {Vinj[i]}-18.75NaOH.csv', index=False)
In [ ]:
```