```
In [1]:
           1 import warnings
           2 import numpy as np
           3 import pandas as pd
           4 from pyIClab import (
                  IonExchanger, Column, Eluent, SwitchingValve, PEEKTubing,
                  SampleLoop, DSM_SimpleEquilibriums, ContaminatedPhreeqcSuppressorBeta,
           7
                  Detector, IonChromatograph, DSM_SEConstrutor
           8
                  )
           9 import seaborn as sns
          10 | import matplotlib.pyplot as plt
          11 from pyIClab.beadedbag import mpl_custom_rcconfig
          12 from IPython.display import clear_output
          13 from scipy.integrate import quad
          14 from scipy.interpolate import interp1d
```

```
In [2]: ▼
            1 test_params = dict(
                   test_name = 'Gradient E',
            2
            3
                   fname = 'as18-20240731.dat',
            4
                   directory = 'db',
            5
                   profile = {
            6
                        'OH-':(
            7
                            (0, 17),
            8
                            (0.01, 17),
            9
                            (5, 38.5),
           10
                            (6, 38.5),
                            (6, 17),
           11
                            (10, 17)
           12
           13
                            ),
           14
                        },
                   length = '25 \text{ cm'},
           15
                   inner_diameter = '4 mm',
           16
           17
                   CO2\_level = .01,
                   sample = {
           18
           19
                        'F-': '0.1 mM',
                        'Cl-': '0.1 mM',
           20
           21
                        'NO2-': '0.2 mM',
           22 #
                          'Br-': '0.2 mM',
                        'NO3-': '0.2 mM',
           23
                        'SO4-2': '0.1 mM',
           24
                        'PO4-3': '0.2 mM',
           25
           26
                        },
           27
                   )
           28
```

```
In [3]:
            1 | from pyIClab import DSM_CEConstrutor
              from pyIClab.engines.equilibriums import find_x_LSSM
            3
            4
              class LocalConstructor(DSM_CEConstrutor):
            5
            6
                  def set_x(self):
            7
            8
                       kmap = self.set_kmap()
            9
                       return find_x_LSSM(kmap, -1)
           10
           11
```


```
In [4]: ▼
           1 | def unit_test(
            3
                  test_name,
                  fname,
            5
                  directory,
            6
                  length,
            7
                  inner_diameter,
            8
                  profile,
            9
                  CO2_level,
           10
                  sample,
           11
                  model_constructor_prompt,
           12
           13
           14
                   # -----
                  sp = IonExchanger.load(fname, directory=directory)
           15
                  column = Column('Column', length=length, ID=inner_diameter)
           16
                  column.pack(sp)
           17
           18
           19
                   # -----
                  eluent = Eluent(name='EG', profile=profile)
           20
                  tb = SampleLoop('PEEK', V='0.21 mL') # hold-up volume
           21
           22
                  sixport = SwitchingValve.SixPort()
           23
                  loop = SampleLoop('Loop', '25 uL')
           24
                  suppressor = ContaminatedPhreeqcSuppressorBeta('Suppressor', 'anion', _CO2_level=CO2_level)
           25
                  detector = Detector('Detector')
           26
                   # -----
           27
           28
                  eluent.assemble(tb)
           29
                  sixport.assemble(0, tb)
           30
                  sixport.assemble(1, column)
           31
                  sixport.assemble([2, 5], loop)
           32
                  column.assemble(suppressor)
           33
                  suppressor.assemble(detector)
           34
                  ic = IonChromatograph('Gradient-Test', ('OH-',), lockon=sixport)
           35
                  commands = '''
           36
                  0.0 min, sixport, inject
           37
           38
                  0.5 min, sixport, load
           39
           40
                  ic.reset_commands(commands)
           41
                  prompt = model_constructor_prompt if model_constructor_prompt != 'DSM_CE_NIC' else LocalConstructor
           42
                  ic.set_ModelConstructor(prompt , column)
           43
           44
                   # -----
                  water = {'Cl-': '1e-9 mM'}
           45
           46
                  ic.inject(water, loop)
           47
                  ic.go(tmax=eluent._tmax)
           48
                  df1 = detector.get_signals(signal_type='conductivity')
           49
                  df1.to_csv(
                      f'''{test_name}-{str(model_constructor_prompt).replace('_', '-')}'''
           50
                       '''-Background.txt''',
           51
                      index=False,
           52
           53
           54
           55
                   # -----
           56
                  ic.inject(sample, loop)
           57
                  ic.go(tmax=eluent._tmax)
           58
                  df2 = detector.get_signals(signal_type='conductivity')
           59
                      f'''{test_name}-{str(model_constructor_prompt).replace('_', '-')}'''
           60
                       '''-Total.txt''
           61
           62
                      index=False,
           63
                      )
           64
                  return df1, df2
           65
              for model_constructor_prompt in ['DSM_SE', 'DSM_CE', 'DSM_CE_NIC']:
In [5]: ▼
            1
                  with warnings.catch_warnings(action='ignore'):
            2
            3
                      df1, df2 = unit_test(
                          model_constructor_prompt=model_constructor_prompt,
            4
```

```
5
               **test_params,
6
           clear_output()
```


```
In [6]:
            1 backgrounds = {}
            2 chroms = {}
              for model in ['DSM-SE', 'DSM-CE', 'DSM-CE-NIC']:
                  test_name = test_params.get('test_name')
            5
                  fname_bg = f'{test_name}-{model}-Background.txt'
            6
                  backgrounds[model] = pd.read_csv(fname_bg)
            7
                  fname_tt = f'{test_name}-{model}-Total.txt'
            8
                   chroms[model] = pd.read_csv(fname_tt)
            9
           10
              df_exp = pd.read_csv(f'{test_name}-exp-Total.txt',
                  sep='\s+',
           11
           12
                  skiprows=43,
                  names=['time', 'step', 'signal'],
           13
           14
                  )[['time', 'signal']]
           15
```

```
In [7]:
            1 sns.set()
            2 plt.rcParams.update(mpl_custom_rcconfig)
            3 fig, ax = plt.subplots()
              for i, model in enumerate(['DSM-SE', 'DSM-CE', 'DSM-CE-NIC']):
            6
                  df = chroms[model]
            7
            8
                  x, y = df['time'], df['signal']
                  ax.plot(x, y, label=model)
           10
           11 | ax.plot('time', 'signal', data=df_exp, label='Experimental', linestyle='--')
           12 | ax.set_xlabel('time, min', fontsize=10, fontweight='bold')
           13 | ax.set_ylabel('conductivity, μS/cm', fontsize=10, fontweight='bold')
           14 ax.set(xlim=(0, max(x)), ylim=(-1, max(y)*1.2))
           15 ax.legend()
           16
          17
             # for i in peaks:
                    tR = x[i]
           18 #
           19
                    signal = y[i]
           20
                    ax.text(tR, signal + .1, f'{tR:.2f}', ha='center', zorder=2)
           21
           22 profile = test_params.get('profile')
           23 ax2 = ax.twinx()
           24 | ax2.plot(*zip(*profile['OH-']), color='grey', linestyle='--', zorder=1)
           25 ax2.grid(visible=False)
           26 | ax2.set_ylabel('[KOH], mM', fontsize=10, fontweight='bold')
```

Out[7]: Text(0, 0.5, '[KOH], mM')



```
In [19]:
             1 df0 = backgrounds['DSM-SE']
             2 df1 = chroms['DSM-SE']
             3 f0 = interp1d(df0['time'], df0['signal'])
             4 f1 = interp1d(df1['time'], df1['signal'])
             5 f = lambda t: f1(t) - f0(t)
             7 \text{ window} = (5.7, 6.7)
             8 t = np.linspace(*window, 10000, endpoint=True)
             9 df = pd.DataFrame(data=dict(time=t, signal=f(t)))
            10 ax = df.plot(x='time', y='signal')
           11 | ax.text(df['time'][df['signal'].argmax()], df['signal'].max(),
                    s=round(df['time'][df['signal'].argmax()], 3),
                   ha='center',
va='bottom',)
            13
            14
           with warnings.catch_warnings(action='ignore'):
                   ax.text(0.1, .5, f'Peak Area: {round(quad(f, *window)[0], 4)}',
           16
            17
                        transform=ax.transAxes,)
            18
```



```
In [ ]: 1
```