

```
In [1]: 1 from pyIClab import (  
2     DSM_CEConstrutor, IonChromatograph, Eluent, SwitchingValve,  
3     Column, IonExchanger, SampleLoop, Detector,  
4 )  
5 from pyIClab import ContaminatedPhreeqcSuppressorBeta as Suppressor  
6 from pyIClab.engines.equilibriums import find_x_LSSM  
7 from pyIClab.engines.models import _total_mix  
8 from pyIClab.beadedbag import mpl_custom_rcconfig  
9 import matplotlib.pyplot as plt  
10 import seaborn as sns  
11 import numpy as np  
12 from IPython.display import clear_output  
13 import warnings  
14  
15 sns.set()  
16 plt.rcParams.update(mpl_custom_rcconfig)  
17
```

```
In [2]: 1 def local_post_distrubute(model, /, *,  
2     mix_n: int,  
3     ):  
4     '''  
5     Just for prettier water dips.  
6     '''  
7  
8     for _ in range(mix_n):  
9         _total_mix(model)
```

```
In [3]: 1 class LocalConstructor(DSM_CEConstrutor):  
2  
3     def set_x(self):  
4  
5         kmap = self.set_kmap()  
6  
7         return find_x_LSSM(kmap, -1)  
8  
9     def set_post_distribute(self):  
10  
11         return local_post_distrubute  
12  
13     def set_post_distribute_params(self):  
14  
15         N = self.set_N()  
16         length = self.host.length.to('cm').magnitude  
17         target_N = round(length / 0.04)  
18         mix_n = round(np.log2(2*N / target_N)) + 1  
19  
20         return {'mix_n': mix_n}
```

```
In [4]: 1 solutions = [{'C1-1': f'{c} mM'} for c in [.1, .5, 2.5, 10, 100]]
```

```
In [5]: 1 ic_collection = []  
2 for i in range(len(solutions)):  
3     eluent = Eluent.HydroxideIsocratic('18.75 mM', name=f'EG{i}')  
4     sp = IonExchanger.load('home_made.dat', directory='db')  
5     column = Column(f'Home_made{i}', length='15 cm', ID='4.6 mm')  
6     column.pack(sp)  
7     sixport = SwitchingValve.SixPort(name=f'Sixport{i}')  
8     suppressor = Suppressor(name=f'Suppressor{i}', kind='anion', _CO2_level=.018)  
9     detector = Detector(name=f'Detector{i}')  
10    loop = SampleLoop(name=f'Loop{i}', V='25 uL')  
11  
12    sixport.assemble(0, eluent)  
13    sixport.assemble([2, 5], loop)  
14    sixport.assemble(1, column)  
15    column.assemble(suppressor)  
16    suppressor.assemble(detector)  
17  
18    ic = IonChromatograph(name=f'IC{i}', competing_ions=('OH[-1]',), lockon=sixport)  
19    ic_collection.append(ic)  
20  
21 ic_collection
```

```
Out[5]: [<IC System "IC0">,  
<IC System "IC1">,  
<IC System "IC2">,  
<IC System "IC3">,  
<IC System "IC4">]
```

```
In [6]: 1 for i, ic in enumerate(ic_collection):
2
3     solution = solutions[i]
4     ic.inject(solution, f'loop{i}')
5     commands = f'''
6         0.0 min, sixport{i}, inject
7     '''
8     ic.reset_commands(commands)
9     ic.set_ModelConstructor(LocalConstructor, f'Home_made{i}')
```

▼

```
10
11 with warnings.catch_warnings(action='ignore'):
12     ic.start(tmax='10 min')
13
14 clear_output()
```

```
In [7]: 1 df_collection = []
2 for i, ic in enumerate(ic_collection):
3
4     (detector,) = tuple(ic.detectors)
5     df_collection.append(detector.get_signals(signal_type='conductivity'))
```

Calculating eluent conductivity on <Detector "Detector0">...: 0%| | 0/5999 [00:00<?, ?it/s]


Calculating eluent conductivity on <Detector "Detector1">...: 0%| | 0/5999 [00:00<?, ?it/s]

Calculating eluent conductivity on <Detector "Detector2">...: 0%| | 0/5999 [00:00<?, ?it/s]

Calculating eluent conductivity on <Detector "Detector3">...: 0%| | 0/5999 [00:00<?, ?it/s]

Calculating eluent conductivity on <Detector "Detector4">...: 0%| | 0/5999 [00:00<?, ?it/s]

```
In [8]: 1 fig, ax = plt.subplots()
2 offsets = [0, 0, 0, 0, 0, 0]
3 for i, df in enumerate(df_collection):
4     ax.plot(df['time'], df['signal']+offsets[i])
5
6 # ax.set(ylim=(0, 1))
```



```
In [9]: 1 from scipy.signal import find_peaks
2 tR = []
3 for df in df_collection:
4     (pk, _), _ = find_peaks(df['signal']-df['signal'][0], height=2)
5     tR.append(df['time'][pk])
6
7 tR
```

Out[9]: [4.952822440813605,  
4.9711600533511175,  
5.056179893297767,  
5.25122540846949,  
6.0197380793597866]

```
In [10]: 1 # for i, (s, df) in enumerate(zip(solutions, df_collection)):
2
3 #     df.to_csv(f'effective charge ratio {s.copy().popitem()}-18.75NaOH.csv', index=False)
```

```
In [ ]: 1
```

